TRAP-6

Cat # 11005 1 flacon à 500 µM

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Cat # 11005 1 vial, 500 µM

Last revision: 2018-06-20

For Research Use Only. Not for use in diagnostic procedures.

INTENDED USE

Thrombin receptor activator peptide-6 (TRAP-6) is a strong agonist that causes platelet aggregation through one of the thrombin receptors (PAR-1 or protease activated receptor-1). When added to citrated platelet-rich plasma (cPRP), TRAP-6 induces a strong and irreversible aggregation which arises from release of platelet granule content and thromboxane A2 synthesis.

CAUTION

For *in vitro* use only.

The disposal of waste materials must be carried out according to current local regulations.

REAGENT PREPARATION

- Reconstitute with distilled water.
For example, reconstitute the vial with exactly 1.0 ml of distilled water to obtain a TRAP-6 solution at 1 mM.
Allow the reconstituted solution to stand for 10 minutes at 18-25°C.
Homogenize by inversion and swirling before use.
- After TRAP-6 solution reconstitution, a working solution will be prepared by dilution in normal saline solution (NaCl 0.9%) or in other isotonic buffer (Tyrode's, PBS...).

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REAGENT STORAGE AND STABILITY

Stored at 2-8°C, the reagent in unopened vial is stable until the expiration date indicated on the vial label.

Once reconstituted, the reagent is stable:

- 8 hours at 18-25°C,
- 2 weeks at 2-8°C,
- 2 months at -20°C.

REAGENTS AND MATERIAL REQUIRED (NOT PROVIDED)

- Appropriate disposables for performing platelet aggregation studies
- Aggregometer
- Normal saline solution (NaCl 0.9%) or Tyrode's buffer or PBS
- Common research laboratory equipment and materials

EXAMPLE OF USE

TRAP-6 induced platelet aggregation test is classically performed on cPRP using a 10 µM final concentration. Other concentration conditions may be used with this reagent.

PERFORMANCES

Within run reproducibility (n=21) in human cPRP using an aggregometer:

Concentration in cPRP	Max. aggregation	CV
10 µM	87 %	4.2 %

BIBLIOGRAPHIE

- Cattaneo M, Cerletti C, Harrison P *et al* (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH, J. Thromb. Haemost. 11, 1183-1189
- Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures Rand with the aim to develop standards in the field. Am. J. Clin. Pathol. 123, 172-183
- Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH, J. Thromb. Haemost. 7, 1029
- Kehrel BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
- Gremmel T, Frelinger III AL, Michelson AD. (2017) Platelet Physiology, Semin. Thromb. Hemost. 42(3):191-204

REFERENCES

- Cattaneo M, Cerletti C, Harrison P *et al* (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH, J. Thromb. Haemost. 11, 1183-1189
- Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures Rand with the aim to develop standards in the field. Am. J. Clin. Pathol. 123, 172-183
- Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH, J. Thromb. Haemost. 7, 1029
- Kehrel BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
- Gremmel T, Frelinger III AL, Michelson AD. (2017) Platelet Physiology, Semin. Thromb. Hemost. 42(3):191-204

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PLATELET RANGE: OTHER AVAILABLE REAGENTS

Denomination	Cat #	Contents
Arachidonic acid	11001	1 vial, 15 mM
ADP	11002	1 vial, 1 mM
Equine collagen	11003	1 vial, 100 µg/ml
Epinephrine	11004	1 vial, 1 mM
Murine collagen	11006	1 vial, 1 mg/ml
Apyrase ADP-High	PY062310	1 vial, 10 U
	PY0623100	1 vial, 100 U
	PY0623200	1 vial, 200 U
	PY0623500	1 vial, 500 U
Apyrase ADP-Premium	PY062710	1 vial, 10 U
	PY0627100	1 vial, 100 U
	PY0627200	1 vial, 200 U
	PY0627500	1 vial, 500 U

15102-03

English

AGRO-BIO
2 allée de la Chavannerie
45240 La Ferté Saint-Aubin
France

+33 (0)2 38 64 83 50
techsupport@agro-bio.com
www.agro-bio.com

Information and/or pictures contained in this document are protected by copyright and other intellectual property rights. © 2017, Agro-Bio, all rights reserved. Agro-Bio's logos and products names are registered trademarks.

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- Cattaneo M, Cerletti C, Harrison P et al (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH, J. Thromb. Haemost. 11, 1183-1189
- Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures Rand with the aim to develop standards in the field. Am. J. Clin. Pathol. 123, 172-183
- Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH, J. Thromb. Haemost. 7, 1029
- Kehrel BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
- Gremmel T, Frelinger III AL, Michelson AD. (2017) Platelet Physiology, Semin. Thromb. Hemost. 42(3):191-204

LITERATUR

- Cattaneo M, Cerletti C, Harrison P et al (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH, J. Thromb. Haemost. 11, 1183-1189
- Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures Rand with the aim to develop standards in the field. Am. J. Clin. Pathol. 123, 172-183
- Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH, J. Thromb. Haemost. 7, 1029
- Kehrel BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
- Gremmel T, Frelinger III AL, Michelson AD. (2017) Platelet Physiology, Semin. Thromb. Hemost. 42(3):191-204

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Les informations et/ou images contenues dans ce document sont protégées par copyright et autres droits de propriétés intellectuelles. © 2017, Agro-Bio, tous droits réservés. Les logos et/ou les noms de produits d'Agro-Bio sont des marques déposées.
Las informaciones y/o imágenes contenidas en este documento están protegidas por copyright y otros derechos de propiedad intelectual. © 2017, Agro-Bio, todos los derechos reservados. Los logos y/o los nombres de productos de Agro-Bio son marcas registradas.
Die in diesem Dokument enthaltenen Informationen und/oder Abbildungen sind durch Copyright und andere Rechte am geistigen Eigentum geschützt. © 2017, Agro-Bio, alle Rechte vorbehalten. Die Logos und/oder Produktnamen von Agro-Bio sind eingetragene Markenzeichen.