APYRASE ADP-HIGH EC 3.6.1.5

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BIBLIOGRAPHIE

- Kalckar HM. (1994) Adenylpyrophosphatase and myokinase. J. Biol. Chem. 153, 355-367
- Kettlun A. *et al.* (1982) Properties of two apyrases from Solanum tuberosum. Phytochemistry 21, 551
- Molnar J. and Lorand L. (1961) Studies on apyrases. Arch. Biochem. Biophys. 93, 353-363
- Abbraccio MP. *et al.* (2009) Purinergic signalling in the nervous system: an overview. Trends Neurosci. Jan 32 (1), 19-29
- Burnstock G. (2009) Purinergic signalling: past, present and future. Braz. J. Med. Biol. Res. Jan 42 (1), 38
- Gachet C. (2011) P2 receptors and platelet function. Purinergic Signal. Sep 7 (3), 293-303
- Burnstock G. & Boeynaems JM. (2014) Purinergic signaling and immune cells. Purinergic Signal. Dec 10 (4), 529-564
- Sévigny J. *et al.* (2015) Purinergic signaling in immune system regulation in health and disease. Mediators of Inflammation. Article ID 106863 Editorial 1-3
- Hechler B. & Gachet C. (2015) Purinergic receptors in thrombosis and inflammation. Arterioscler. Thromb. Vasc. Biol. Sep 10, 1-9
- Ferrari D. *et al.* (2015) Purinergic signaling in scarring. FASEB J. Sep (2) Published ahead
- Burnstock G. & Pelleg A. (2015) Cardiac purinergic signaling in health and disease. Purinergic Signal. Mar 11 (1), 1-46
- Cazenave JP. *et al.* (2004) Preparation of washed platelet suspensions from human and rodent blood. Methods Mol. Biol. 272, 13-28
- Baurand A. *et al.* (2000) Desensitization of the platelet aggregation response to ADP, differential down-regulation of the P2Y₁ and P2cy₂ receptors. Thromb. Haemost. 85, 303-308
- Vial C. *et al.* (1997) Presence of P2X₁ purinoceptors in human platelets and megakaryoblastic cell lines. Thromb. Haemost. 78 (6), 1500-1504
- Mahaut Smith MP. *et al.* (2004) Emerging roles for P2X₁ receptors in platelet activation. Platelets 15 (3), 131-144
- Kauffmanstein G. *et al.* (2014) Nucleoside triphosphate diphosphohydrolase-1 ectonucleotidase is required for normal vas deferens contraction and male fertility through maintaining P2X₁ receptor function. J. Biol. Chem. 289 (41), 28629-28639


APYRASE ADP-HIGH EC 3.6.1.5

REF	Presentation		
PY062310	1 vial	10 U	lyophilized
PY0623100	1 vial	100 U	lyophilized
PY0623200	1 vial	200 U	lyophilized
PY0623500	1 vial	500 U	lyophilized
On request	1 vial	1000 U	frozen
On request	1 vial	2000 U	frozen
On request	1 vial	5000 U	frozen

FOR RESEARCH USE ONLY.
NOT FOR USE IN DIAGNOSTIC PROCEDURES.

PRODUCT DESCRIPTION

Apyrase is an enzyme extracted from a selected variety of *Solanum tuberosum* able to catalyze the sequential hydrolysis of ATP to ADP and ADP to AMP, while releasing inorganic phosphate molecules (Pi)¹⁻². It exists in more than one molecular form and at least two isoenzymes have been identified, which differ in their relative rates of hydrolysis of ADP and ATP²⁻³. **APYRASE ADP-HIGH** has a high ADPase/ATPase ratio.

PRESENTATION

1 vial of purified and lyophilized apyrase with an ADPase activity of 20 to 200 U/mg. The 1000 U, 2000 U and 5000 U references are available solely in the frozen format.

EXCIPIENT

Sodium chloride

RECONSTITUTION

Reconstitute the vial with distilled water.

For example, reconstitute the 100 U vial (*ref.* PY0623100) with exactly 1 mL of distilled water to obtain a solution at 100 U/mL. Allow the reconstituted solution to stand for 15 minutes at 2-8°C. Homogenize by inversion and swirling before use.

REAGENT STORAGE AND STABILITY

Storage: before reconstitution, store at 2-8°C, away from light, until the expiration date indicated on the vial or on the Quality Control flyer.

Stability: after reconstitution, store at 2-8°C. After use, recap and quickly place the vial at 2-8°C. Avoid repeated freezing/thawing cycles.

USE EXAMPLE

APYRASE ADP-HIGH is used in all applications mainly targeting the purinergic P2 receptors including the metabotropic receptors P2Y₁, P2Y₁₂ and P2Y₁₃ in research domains such as immunology, neurology, oncology and hematology⁴⁻¹¹.

APYRASE ADP-HIGH is especially suitable for the preparation of isolated platelets. Adding **APYRASE ADP-HIGH** (e.g. 0.02 U/mL) to a platelet suspension at 100,000 to 300,000 platelets/ μ L preserves platelets sensitivity to ADP¹². In the absence of **APYRASE ADP-HIGH**, the platelets become refractory to ADP and therefore unable to change shape or aggregate. This is due in particular to the release of nucleotides during the washing of platelets, which causes desensitization of ADP receptors P2Y₁₃¹³. **APYRASE ADP-HIGH** (0.2 to 2.0 U/mL) can also be used to **reversibly** inhibit platelet aggregation induced by 5.0 μ M ADP through P2Y₁₂ and P2Y₁ receptors or to target P2X₁ receptors¹⁴⁻¹⁶.

PERFORMANCE CHARACTERISTICS

Since the ADPase/ATPase ratio could vary from one lot to another, it is clearly indicated for each lot on the Quality Control flyer.

An apyrase unit corresponds to the release of 1.0 μ mole of inorganic phosphate from ADP or ATP per minute at pH 6.5 at 37°C.

REFERENCES

- Kalckar HM. (1994) Adenylpyrophosphatase and myokinase. J. Biol. Chem. 153, 355-367
- Kettlun A. *et al.* (1982) Properties of two apyrases from Solanum tuberosum. Phytochemistry 21, 551
- Molnar J. and Lorand L. (1961) Studies on apyrases. Arch. Biochem. Biophys. 93, 353-363
- Abbraccio MP. *et al.* (2009) Purinergic signalling in the nervous system: an overview. Trends Neurosci. Jan 32 (1), 19-29
- Burnstock G. (2009) Purinergic signalling: past, present and future. Braz. J. Med. Biol. Res. Jan 42 (1), 38
- Gachet C. (2011) P2 receptors and platelet function. Purinergic Signal. Sep 7 (3), 293-303
- Burnstock G. & Boeynaems JM. (2014) Purinergic signaling and immune cells. Purinergic Signal. Dec 10 (4), 529-564
- Sévigny J. *et al.* (2015) Purinergic signaling in immune system regulation in health and disease. Mediators of Inflammation. Article ID 106863 Editorial 1-3
- Hechler B. & Gachet C. (2015) Purinergic receptors in thrombosis and inflammation. Arterioscler. Thromb. Vasc. Biol. Sep 10, 1-9
- Ferrari D. *et al.* (2015) Purinergic signaling in scarring. FASEB J. Sep (2) Published ahead
- Burnstock G. & Pelleg A. (2015) Cardiac purinergic signaling in health and disease. Purinergic Signal. Mar 11 (1), 1-46
- Cazenave JP. *et al.* (2004) Preparation of washed platelet suspensions from human and rodent blood. Methods Mol. Biol. 272, 13-28
- Baurand A. *et al.* (2000) Desensitization of the platelet aggregation response to ADP, differential down-regulation of the P2Y₁ and P2cy₂ receptors. Thromb. Haemost. 85, 303-308
- Vial C. *et al.* (1997) Presence of P2X₁ purinoceptors in human platelets and megakaryoblastic cell lines. Thromb. Haemost. 78 (6), 1500-1504
- Mahaut Smith MP. *et al.* (2004) Emerging roles for P2X₁ receptors in platelet activation. Platelets 15 (3), 131-144
- Kauffmanstein G. *et al.* (2014) Nucleoside triphosphate diphosphohydrolase-1 ectonucleotidase is required for normal vas deferens contraction and male fertility through maintaining P2X₁ receptor function. J. Biol. Chem. 289 (41), 28629-28639

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 **AGRO-BIO**

+33 (0)2 38 64 83 50
techsupport@agro-bio.com
www.agro-bio.com

Information and/or pictures contained in this document are protected by copyright and other intellectual property rights. © 2015, Agro-Bio, all rights reserved. Agro-Bio's logos and products names are registered trademarks.

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Tous droits réservés. Les logos et/ou les images contenues dans ce document sont protégées par copyright et autres droits de propriétés intellectuelles. © 2015, Agro-Bio, tous droits réservés. Les logos et/ou les noms de produits d'Agro-Bio sont des marques déposées.