ADP

Cat # 11002

1 flacon à 1 mM

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Cat # 11002

1 vial, 1 mM

For Research Use Only. Not for use in diagnostic procedures.

INTENDED USE

Adenosine diphosphate (ADP) is a weak agonist that causes platelet aggregation through P2Y₁ and P2Y₁₂ receptors activation. When added to citrated platelet-rich plasma (cPRP), ADP induces a two wave aggregation, where the secondary wave arises from release of platelet granule content and thromboxane A₂ synthesis.

CAUTION

For *in vitro* use only.
The disposal of waste materials must be carried out according to current local regulations.

REAGENT PREPARATION

- Reconstitute with distilled water.
For example, reconstitute the vial with exactly 0.5 mL of distilled water to obtain an ADP solution at 1 mM.
Allow the reconstituted solution to stand for 10 minutes at 18-25°C.
Homogenize by inversion and swirling before use.
- After ADP solution reconstitution, a working solution will be prepared by dilution in normal saline solution (NaCl 0.9%) or in other isotonic buffer (Tyrode's, PBS...).

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REAGENT STORAGE AND STABILITY

Stored at 2-8°C, the reagent in unopened vial is stable until the expiration date indicated on the vial label.

Once reconstituted, the reagent is stable:

- 2 weeks at 2-8°C,
- 8 hours at 18-25°C.

REAGENTS AND MATERIAL REQUIRED (NOT PROVIDED)

- Appropriate disposables for performing platelet aggregation studies
- Aggregometer
- Normal saline solution (NaCl 0.9%) or Tyrode's buffer or PBS
- Common research laboratory equipment and materials

EXAMPLE OF USE

ADP induced platelet aggregation test is classically performed on cPRP using a 2 µM final concentration. Other concentration conditions may be used with this reagent.


Figure 1. Platelet response to different concentrations of ADP on cPRP using a commercial aggregometer (brand A)

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
Figure 2. Platelet response to different concentrations of ADP on cPRP using a commercial aggregometer (brand B)


Figure 3. Platelet response to different concentrations of ADP on cPRP using a commercial aggregometer (brand C)

PERFORMANCES

Repeatability at a 2 µM final concentration using three aggregometers of brands A, B and C, respectively.

Aggregometer	n	Max. aggregation (%)	CV (%)
A	10	90	12,0
B	10	85	9,2
C	10	77	5,8

BIBLIOGRAPHIE

- Cattaneo M, Cerletti C, Harrison P et al (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH, J. Thromb. Haemost. 11, 1183-1189
- Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures Rand with the aim to develop standards in the field. Am. J. Clin. Pathol. 123, 172-183
- Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH. J. Thromb. Haemost. 7, 1029
- Köhnel BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
- Grammel T, Frelinger III AL, Michelson AD. (2016) Platelet Physiology. Semin. Thromb. Hemost. 42(3):191-204

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REFERENCES

- Cattaneo M, Cerletti C, Harrison P et al (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH, J. Thromb. Haemost. 11, 1183-1189
- Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures Rand with the aim to develop standards in the field. Am. J. Clin. Pathol. 123, 172-183
- Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH. J. Thromb. Haemost. 7, 1029
- Köhnel BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
- Grammel T, Frelinger III AL, Michelson AD. (2016) Platelet Physiology. Semin. Thromb. Hemost. 42(3):191-204

PLATELET RANGE: OTHER AVAILABLE REAGENTS

Denomination	Cat #	Contents
Murine collagen	11006	1 vial, 1 mg/mL
Equine collagen	11003	1 vial, 1 mg/mL
Epinephrine	11004	1 vial, 1 mM
Arachidonic acid	11001	1 vial, 50 mM
TRAP-6	11005	1 vial, 1 mM
Apyrase ADP-High	PY062310	1 vial, 10 U
	PY0623100	1 vial, 100 U
	PY0623200	1 vial, 200 U
	PY0623500	1 vial, 500 U
Apyrase ADP-Premium	PY062710	1 vial, 10 U
	PY0627100	1 vial, 100 U
	PY0627200	1 vial, 200 U
	PY0627500	1 vial, 500 U

15086-02 August 2016

English

AGRO-BIO

2 allée de la Chavannerie
45240 La Ferté Saint-Aubin
France

+33 (0)2 38 64 83 50
techsupport@agro-bio.com
www.agro-bio.com

Information and/or pictures contained in this document are protected by copyright and other intellectual property rights. © 2016, Agro-Bio, all rights reserved. Agro-Bio's logos and products names are registered trademarks.

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1. Cattaneo M, Cerletti C, Harrison P et al (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH, J. Thromb. Haemost. 11, 1183-1189
2. Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures and the aim to develop standards in the field. Am. J. Clin. Pathol. 123, 172-183
3. Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH. J. Thromb. Haemost. 7, 1029
4. Kehl BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
5. Grammel T, Freilinger III AL, Michelson AD. (2016) Platelet Physiology, Semin. Thromb. Hemost. 42(3):191-204

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Les informations et/ou images contenues dans ce document sont protégées par copyright et autres droits de propriétés intellectuelles. © 2016, Agro-Bio, tous droits réservés. Les logos et/ou les noms de produits d'Agro-Bio sont des marques déposées.
Las informaciones y/o imágenes contenidas en este documento están protegidas por copyright y otros derechos de propiedad intelectual. © 2016, Agro-Bio, todos los derechos reservados. Los logos y/o los nombres de productos de Agro-Bio son marcas registradas.
Die in diesem Dokument enthaltenen Informationen und/oder Abbildungen sind durch Copyright und andere Rechte am geistigen Eigentum geschützt. © 2016, Agro-Bio, alle Rechte vorbehalten. Die Logos und/oder Produktnamen von Agro-Bio sind eingetragene Markenzeichen.