## Epinéphrine

Cat # 11004

1 flacon à 1 mM

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Cat # 11004

1 vial, 1 mM

For Research Use Only. Not for use in diagnostic procedures.

### INTENDED USE

Epinephrine is a weak agonist that causes platelet aggregation through its effect on receptor  $\alpha_2A$ . When added to citrated platelet-rich plasma (cPRP), epinephrine induces a two wave aggregation, where the secondary wave arises from release of platelet granule content and thromboxane A2 synthesis.

### CAUTION

For *in vitro* use only.

The disposal of waste materials must be carried out according to current local regulations.

### REAGENT PREPARATION

- 1- Reconstitute with distilled water.

For example, reconstitute the vial with exactly 0.5 mL of distilled water to obtain an epinephrine solution at 1 mM.

Allow the reconstituted solution to stand for 10 minutes at 18-25°C. Homogenize by inversion and swirling before use.

- 2- After epinephrine solution reconstitution, a working solution will be prepared by dilution in normal saline solution (NaCl 0.9%) or in other isotonic buffer (Tyrode's, PBS...).

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### REAGENT STORAGE AND STABILITY

Stored at 2-8°C, the reagent in unopened vial is stable until the expiration date indicated on the vial label.

Once reconstituted, the reagent is stable:

- 2 weeks at 2-8°C,
- 8 hours at 18-25°C,

### REAGENTS AND MATERIAL REQUIRED (NOT PROVIDED)

- Appropriate disposables for performing platelet aggregation studies
- Aggregometer
- Normal saline solution (NaCl 0.9%) or Tyrode's buffer or PBS
- Common research laboratory equipment and materials

### EXAMPLE OF USE

Epinephrine induced platelet aggregation test is classically performed on cPRP using a 5  $\mu$ M final concentration. Other concentration conditions may be used with this reagent.


Figure 1. Platelet response to 5  $\mu$ M of epinephrine on cPRP using a commercial aggregometer (brand A)

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
Figure 2. Platelet response to 5  $\mu$ M of epinephrine on cPRP using a commercial aggregometer (brand B)


Figure 3. Platelet response to 5  $\mu$ M of epinephrine on cPRP using a commercial aggregometer (brand C)

### PERFORMANCES

Repeatability at a 5  $\mu$ M final concentration using three aggregometers of brands A, B and C, respectively.

| Aggregometer | n  | Max. aggregation (%) | CV (%) |
|--------------|----|----------------------|--------|
| A            | 10 | 101                  | 10.4   |
| B            | 10 | 81                   | 8.9    |
| C            | 10 | 78                   | 7.1    |

### BIBLIOGRAPHIE

1. Cattaneo M, Cerletti C, Harrison P et al (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH. J. Thromb. Haemost. 11, 1183-1189
2. Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures Rand with the aim to develop standards in the field. Am. J. Clin. Pathol. 123, 172-183
3. Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH. J. Thromb. Haemost. 7, 1029
4. Kehl BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
5. Grammel T, Freilinger III AL, Michelson AD. (2016) Platelet Physiology. Semin. Thromb. Hemost. 42(3):191-204

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### REFERENCES

1. Cattaneo M, Cerletti C, Harrison P et al (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH. J. Thromb. Haemost. 11, 1183-1189
2. Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures Rand with the aim to develop standards in the field. Am. J. Clin. Pathol. 123, 172-183
3. Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH. J. Thromb. Haemost. 7, 1029
4. Kehl BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
5. Grammel T, Freilinger III AL, Michelson AD. (2016) Platelet Physiology. Semin. Thromb. Hemost. 42(3):191-204

### PLATELET RANGE: OTHER AVAILABLE REAGENTS

| Denomination        | Cat #     | Contents        |
|---------------------|-----------|-----------------|
| Murine collagen     | 11006     | 1 vial, 1 mg/mL |
| ADP                 | 11002     | 1 vial, 1 mM    |
| Equine collagen     | 11003     | 1 vial, 1 mg/mL |
| Arachidonic acid    | 11001     | 1 vial, 50 mM   |
| TRAP-6              | 11005     | 1 vial, 1 mM    |
| Apyrase ADP-High    | PY062310  | 1 vial, 10 U    |
|                     | PY0623100 | 1 vial, 100 U   |
|                     | PY0623200 | 1 vial, 200 U   |
|                     | PY0623500 | 1 vial, 500 U   |
| Apyrase ADP-Premium | PY062710  | 1 vial, 10 U    |
|                     | PY0627100 | 1 vial, 100 U   |
|                     | PY0627200 | 1 vial, 200 U   |
|                     | PY0627500 | 1 vial, 500 U   |

15115-02 August 2016

English

AGRO-BIO

2 allée de la Chavannerie  
45240 La Ferté Saint-Aubin  
France

+33 (0)2 38 64 83 50  
techsupport@agro-bio.com  
www.agro-bio.com

Information and/or pictures contained in this document are protected by copyright and other intellectual property rights. © 2016, Agro-Bio, all rights reserved. Agro-Bio's logos and products names are registered trademarks.

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1. Cattaneo M, Cerletti C, Harrison P et al (2013) Recommendations for the standardization of light transmission aggregometry: a consensus of the working party from the platelet physiology subcommittee of SSC/ISTH. J. Thromb. Haemost. 11, 1183-1189
2. Zhou L, Schmaier AH (2005) Platelet aggregation testing in platelet-rich plasma: description of procedures and standards in the field. Am. J. Clin. Pathol. 123, 172-183
3. Cattaneo M, Hayward CP, Moffat KA, Pugliano MT, Liu Y, Michelson AD. (2009) Results of a worldwide survey on the assessment of platelet function by light transmission aggregometry: a report from the platelet physiology subcommittee of the SSC of the ISTH. J. Thromb. Haemost. 7, 1029
4. Kehrel BE, Brode MF. (2013) State of the art in platelet function testing. Transfus. Med. Hemother. 40, 73-86
5. Grammel T, Freilinger III AL, Michelson AD. (2016) Platelet Physiology, Semin. Thromb. Hemost. 42(3):191-204

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Las informaciones y/o imágenes contenidas en este documento están protegidas por copyright y otros derechos de propiedad intelectual. © 2016, Agro-Bio, todos los derechos reservados. Los logos y/o los nombres de productos de Agro-Bio son marcas registradas.
Die in diesem Dokument enthaltenen Informationen und/oder Abbildungen sind durch Copyright und andere Rechte am geistigen Eigentum geschützt. © 2016, Agro-Bio, alle Rechte vorbehalten. Die Logos und/oder Produktnamen von Agro-Bio sind eingetragene Markenzeichen.